# 12-month results overview



#### Conclusions

- At 30 days clinical results show **Major Adverse Event (MAE) composite** of **0.0** % for the Passeo-18 Lux Drug-Coated Balloon (DCB) vs. 8.3 % compared to the control Percutaneus Transluminal Angioplasty (PTA) balloon
- At 6 months angiographic follow-up, Passeo-18 Lux demonstrated a Target Lesion Primary Patency (TLP) of 82.9 % vs. 73.9 % compared to the control PTA balloon
- At 6 months, **59** % of patients **improved in Rutherford Classification** in the DCB group vs. 47 % in the control group. **Improvement of Rutherford Class 5** patients at 6 months was significant in the DCB group (**p = 0.002\***)
- The Passeo-18 Lux is safe in infrapopliteal lesions demonstrated in a low amputation rate and no additional amputations beyond 180 days

(p < 0.05 = significant)

### Study design

#### Design

Prospective, multi-center, randomized controlled, First-In-Man study to assess the safety and performance of the Passeo-18 Lux Drug-Coated Balloon vs. the uncoated Passeo-18 balloon catheter in patients with stenosis and occlusion of the infrapopliteal arteries.

#### Principal investigator


Prof. T. Zeller, Universitäts-Herzzentrum Freiburg, Bad Krozingen, Germany.

#### Primary endpoint

30-day Major Adverse Event<sup>1</sup> rate. 6-month Target Lesion Primary Patency measured by Quantitative Vascular Angiography<sup>2</sup> (QVA).

#### Secondary endpoints

6-month change in Rutherford Class 12-month Major Amputation rate



## Key baseline demographics

Treated length (mm)<sup>2</sup>
Diameter stenosis (%)<sup>2</sup>

Rutherford Class 5 (n/%)

	DCB		PTA	
	Lesion N = 50		Lesion N = 54	
	Mean ± SD	(Min - Max)	Mean ± SD	(Min - Max)
	113.1 ± 88.1	(24.0 to 350.6)	115.0 ± 86.9	(39.2 to 295.0)
	72.5 ± 25.4	(31.0 to 100.0)	72.1 ± 23.2	(30.0 to 100.0)
	<b>DCB</b> Patients N = 36		<b>PTA</b> Patients N = 36	
	26/72 2		24/72 2	


<sup>&</sup>lt;sup>1</sup> MAE = all cause death, major amputation of target extremity, TLR, TVR, target lesion thrombosis, adjudicated by an independent Clinical Events Committee



 $<sup>^{\</sup>rm 2}$  Assessed by an independent Core Laboratory


# 30-day Major Adverse Events (MAE)

(adjudicated by an independent Cinical Events Committee)




## 6-month Target Lesion Primary Patency (TLP)


(assessed by an independent Core Laboratory)



# 6-month change in Rutherford Classification (RC)



#### PTA at 6 months: Improvement: 47 %



#### ■ 12-month Major Amputations

(adjudicated by an independent Clinical Events Committee)

