



BIOSOLVE-II

BIOSOLVE-II 6-month follow-up

Conclusions

- In-segment Late Lumen Loss of 0.27 mm.
- BIOSOLVE-II demonstrates excellent clinical results with a low TLF rate of 3.3 %.
- No definite nor probable scaffold thrombosis occurred up to 6 months.
- BIOSOLVE-II shows a significant improvement in In-segment Late Lumen Loss compared to BIOSOLVE-I and PROGRESS.

Study design

Prospective, multi-center, First-in-Man, patients with single de novo coronary artery lesions in up to two coronary arteries and a Reference Vessel Diameter (RVD) between 2.2-3.7 mm and \leq 21 mm length

Principal investigator

- Prof. Michael Haude, Neuss, Germany


Endpoints

Primary endpoint

- In-segment Late Lumen Loss (LLL) at 6-month

Secondary endpoints (selected)


- Target Lesion Failure (TLF*) defined as a composite of
 - Cardiac death
 - Target vessel myocardial infarction (MI)
 - Clinically driven Target Lesion Revascularization (TLR)
 - Coronary Artery Bypass Surgery (CABG)




* Society of Cardiovascular Angiography and Interventions (SCAI) definition

Primary endpoint results

6-month follow-up



TLF components




Scaffold thrombosis (definite or probable)

	%	95 % confidence interval
Scaffold thrombosis	0.0	0.0-3.1

There was no definite nor probable scaffold thrombosis up to 6 months.

In-segment Late Lumen Loss comparison BIOSOLVE-II, BIOSOLVE-I and PROGRESS

BIOSOLVE-II shows a significant decrease in in-segment LLL compared to BIOSOLVE-I (6-month follow-up) and PROGRESS (4-month follow-up).



Not currently available in the United States.

BIOTRONIK AG
 Ackerstrasse 6
 8180 Bülach · Switzerland
 Tel +41 (0)44 8645111
 Fax +41 (0)44 8645005
 info.vi@biotronik.com
 www.biotronik.com

© BIOTRONIK AG
 All rights reserved. Specifications
 are subject to modification,
 revision and improvement.

BIOTRONIK
 excellence for life