

# 9-month primary endpoint results of the Orient trial, RCT Orsiro vs. Resolute Integrity

## Conclusions

- In this 375 patient randomized trial with an angiographic primary endpoint of late lumen loss (LLL) at 9 months, Orsiro demonstrated non-inferiority to Resolute Integrity. (Orsiro  $0.10 \pm 0.35$  mm, Resolute Integrity  $0.16 \pm 0.39$  mm, p for non-inferiority  $<0.001$ )
- Orsiro, with its ultrathin struts and bioabsorbable polymer, additionally showed numerically better results for the secondary clinical endpoint, target lesion failure (TLF) out to 12 months. (Orsiro 2.4 %, Resolute Integrity 3.3 %,  $p = 0.0623$ )
- These highly encouraging results reconfirm those of previous Orsiro trials and adds to the solid foundation of clinical evidence that supports the use of Orsiro across a broad range of indications

## Study design

A prospective, multi-center, randomized, controlled trial comparing the Orsiro hybrid DES to Resolute Integrity

### Principal Investigator

- Prof. Tae-Jin Youn, Seoul National University Bundang Hospital, South Korea


## Endpoints

### Primary endpoint

- In-stent LLL at 9 months


### Secondary endpoints (selected)

- TLF defined as composite of cardiac death, TLR and target vessel-related MI
- All-cause and cardiac deaths
- Clinically-driven TLR
- MI (target or non-target vessel-related)
- Definite or probable stent thrombosis (ST)




## Primary endpoint

### In-stent Late Lumen Loss



## Secondary endpoint

### Target Lesion Failure



Source: Kang S. EuroPCR 2016. Oral presentation.