# 3Flow

---

Aspiration Catheter

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Aspirationskateter

Aspiratiekatheter

Imukatetri

Catheter d'aspiration

Aspirationskatheter

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Aspirationskateter

Aspirasyon Kateteri

**Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.**

---

<b>English</b>		
Instructions for Use	page	3
<b>中文</b>		
使用说明	页码	4
<b>Česky</b>		
Návod k použití	strana	5
<b>Dansk</b>		
Brugsanvisning	side	6
<b>Nederlands</b>		
Gebruiksaanwijzing	bladzijde	7
<b>Suomi</b>		
Käyttöohjeet	sivu	8
<b>Français</b>		
Mode d'emploi	page	9
<b>Deutsch</b>		
Gebrauchsanweisung	Seite	10
<b>Ελληνικά</b>		
Οδηγίες χρήσης	σελίδα	12
<b>Magyar</b>		
Használati utasítás	oldal	13
<b>Italiano</b>		
Istruzioni per l'uso	pagina	14
<b>Norsk</b>		
Bruksanvisning	side	16
<b>Polski</b>		
Instrukcja użycia	strona	17
<b>Português</b>		
Instruções de utilização	página	18
<b>Русский</b>		
Инструкция по применению	страница	19
<b>Español</b>		
Instrucciones de uso	página	21
<b>Svenska</b>		
Bruksanvisning	sid	22
<b>Türkçe</b>		
Kullanma Talimatı	sayfa	23

---

## English

CAREFULLY READ ALL INSTRUCTIONS PRIOR TO USE. OBSERVE ALL WARNINGS AND PRECAUTIONS NOTED THROUGHOUT THESE INSTRUCTIONS. FAILURE TO DO SO MAY RESULT IN COMPLICATIONS.

Sterilized with ethylene oxide gas. Non pyrogenic. For one procedure only. Do not re-sterilize. Do not use opened or damaged packages. Destroy product after use. Store in a dry place below 40° C, keep away from light.

## Description

The 3Flow Aspiration Catheter is a single-use design, dual lumen catheter. It has a distal radiopaque tip marker and proximal luer-lock port, the proximal luer-lock port is for connection of the Aspiration Line and Aspiration Syringe.

## Indications and intended use

The 3Flow Aspiration Catheter is indicated for use in the central and peripheral circulatory system, including saphenous vein grafts to:

- Contain and aspirate embolic material (thrombus/debris) while performing percutaneous transluminal coronary angioplasty, percutaneous transluminal angioplasty and/or stenting procedures.

## Warnings

- **Prior to use the packaging and product should be inspected for signs of damage. Never use a damaged product or a product from a damaged package.**
- **The 3Flow Aspiration Catheter should be handled carefully. Prior to use inspect the 3Flow Aspiration Catheter carefully for bends, kinks, or other damage. Do not use a damaged 3Flow Aspiration Catheter.**
- **The 3Flow Aspiration Catheter is designed and intended for single use. Check that the sterile pouch is not damaged before use. Do not re-sterilize or re-use, destroy product after use. Re-use & re-sterilization may compromise the device performance.**

## Precautions

- The 3Flow Aspiration Catheter and accessories should be used in conjunction with fluoroscopic guidance and proper anticoagulation agents.
- This device is designed and intended for single patient use only. DO NOT re-sterilize and/or re-use it.
- As in any elective coronary intervention, it is recommended that the patient has a mean systolic blood pressure greater than or equal to 90 mm Hg in concomitant of IV pressors or intra-aortic balloon pump augmentation.

As with most percutaneous interventions, other potential adverse events include: death, myocardial infarction, coronary or bypass graft thrombosis or occlusion, myocardial ischemia, stroke/CVA, emergent or non-emergent bypass graft surgery, emboli (air, tissue or thrombotic), dissection, arterial perforation, arterial rupture, ventricular fibrillation, hemorrhage, hypotension, pseudoaneurysm at access site, arteriovenous fistula, infection at access site, other bleeding complications at access site.

## Materials required for use with the 3Flow aspiration and infusion system

- Minimum 6F or 7F arterial or venous introducer sheath and minimum 6F or 7F femoral guiding catheter in the appropriate configuration to cannulate the vessel (preferably with side holes if ostial narrowing is present or the guide catheter is occlusive).
- Push pull or rotating hemostatic valve
- Heparinized normal saline.

## Recommendation

- If the use of a second guidewire in addition to the 3Flow Aspiration Catheter is required, please note that an increased size of guiding catheter shall be used in order to avoid friction issues.
- Example: for a 6F 3Flow Aspiration Catheter, use a 7F guiding catheter.

## Directions for use

### Preparation of the 3Flow Aspiration Catheter

The 3Flow Aspiration Catheter is supplied with an Aspiration Line, one locking Aspiration Syringe and 2 filters.

01. Remove the 3Flow Aspiration Catheter and accessories from the package.
02. Fill the Aspiration Syringe with approximately 5-10 ml of heparinized saline and attach the Aspiration Line and Aspiration Syringe to the 3Flow Aspiration Catheter.
03. Open the stopcock on the Aspiration Line and flush the entire length of the 3Flow Aspiration Catheter using all of the heparinized saline contained in the Aspiration Syringe. Close the stopcock.
04. Verify that the stopcock on the Aspiration Line is in the closed position. Retract the plunger of the Aspiration Syringe and pull until it locks at the fully extended position. The 3Flow Aspiration Catheter is completely prepared and is ready for use.
05. Perform aspiration using the 3Flow Aspiration Catheter:
  - a. Load and advance the prepared 3Flow Aspiration Catheter over the guidewire to the tip of the guiding catheter.
  - b. Under fluoroscopy advance the 3Flow Aspiration Catheter and position the distal tip marker proximal above the embolic particles. Stop advancement of the 3Flow Aspiration Catheter if any resistance is encountered.
  - c. Begin aspiration by opening the stopcock on the Aspiration Line. Slowly retract the 3Flow Aspiration Catheter towards the guiding catheter. Blood will enter the Aspiration Syringe until the entire vacuum is gone (or the Aspiration Syringe is filled).

### Notes:

- Should blood not begin filling the syringe within 5 seconds, check the guiding catheter tip placement. Unseat the guiding catheter if necessary to resume flow.
- If no blood is aspirated as a result of unseating the guiding catheter, turn the stopcock off and remove the 3Flow Aspiration Catheter from the patient, either flush the aspiration lumen or use a new 3Flow Aspiration Catheter. DO NOT flush the system while the catheter is still inside the patient vasculature.
- During performance testing, the 3Flow Aspiration Catheter demonstrated the ability to evacuate fluid and debris at a minimum rate of 1 ml/second.
  - d. After completing the aspiration process turn the Aspiration Line stopcock to close off the Aspiration Syringe.

06. Remove the 3Flow Aspiration Catheter

**NOTE: Especially for peripheral intervention (or non-coronary): the guidewire must be withdrawn simultaneously with the 3Flow Aspiration Catheter to avoid any loop creation.**

- a. Withdraw the 3Flow Catheter
- b. Slowly retract and remove the 3Flow Aspiration Catheter. If necessary, loosen the Tuohy Borst of the hemostatic valve to allow easy withdrawal of the distal shaft.

**Note:** Remove the Aspiration Syringe and re-flush the aspiration lumen and wire lumen of the 3Flow Aspiration Catheter with heparinized saline before each reuse. Empty the Aspiration Syringe, reattach to the Aspiration Line, close the stopcock, and retract the plunger to the fully extended lock position.

07. Remove the catheters and follow standard hospital practice for management of the insertion site.

## How supplied

The 3Flow Catheter is supplied with an Aspiration Line, one locking Aspiration Syringe and 2 filters; it is supplied sterile and non-pyrogenic in a tray in its own pouch. The 3Flow Aspiration Catheter including all of its components are intended to be used for one procedure only.

## Warranty/Liability

The product and each component of its system (hereinafter "the product") have been designed, manufactured, tested and packaged with all reasonable care. However, the manufacturer has no control over the conditions under which the product is used and a disturbance of the intended function of the product may occur for various reasons. In this respect, the warnings in this instructions for use are expressly to be considered as an integral part of this disclaimer and provide more detailed information. For this reason, the manufacturer disclaims all warranties, expressed or implied regarding the product, including but not limited to, any warranty of merchantability or fitness for a particular purpose of the product. Product descriptions or user guidelines in publications do not constitute any expressed representation or any expressed or implied warranty. The manufacturer is not liable for any direct, incidental or consequential damages or medical expenses caused by any use, defect, failure or malfunction of the product whether the claim is based on contract, warranty, tort or otherwise. This does not apply in the case of intention or in the case of gross negligence of legal representatives or executive staff of the manufacturer. In commercial transactions relating to merchants, the liability is limited to the compensation of typical damages; compensation for any untypical or incidental damage is excluded. These limitations of liability and warranty are not intended to contravene any mandatory provisions of law applicable in the respective country. If any clause of the disclaimer is considered by a competent court to be invalid or to be in conflict with the applicable law, the remaining part of it shall not be affected and remain in full force and effect. The invalid clause shall be substituted by a valid clause which best reflects the manufacturer's legitimate interest in limiting its liability or warranty without infringing any mandatory provisions of applicable law. No person has any authority to bind the manufacturer to any warranty or liability regarding the product.

## Symbol legend

	Sterilized using ethylene oxide		Do not use if package is damaged
	Do not reuse		Batch code
	Caution, consult accompanying documents		Catalogue number
	Keep dry		Use by
	Keep away from sunlight		Instruction for use
	Temperature limitation		Manufacturer
	Do not re-sterilize		Diameter

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All rights reserved. Specifications subject to
BIOTRONIK, Inc.
BIOTRONIK AG