

ProMRI® System Technical Manual



CAUTION

Federal (U.S.A.) law restricts this device to sale by, or on the order of,
a physician (or properly licensed practitioner).

© 2016 BIOTRONIK, Inc., all rights reserved.

Contents

1. Basic Information	1
1.1 About this Manual	1
1.1.1 Subject of this manual	1
1.1.2 What this manual doesn't include	1
1.2 Target Group	1
1.2.1 Cooperation between professionals from two areas of expertise	1
1.2.2 Knowledge required by the cardiology specialist	1
1.2.3 Knowledge required by the radiology specialist	1
1.3 Active Device and Lead	2
1.3.1 Patient selection, MRI indication	2
1.3.2 Intended use	2
1.3.3 Residual risk	2
1.4 ProMRI® System	2
2. Safety Warnings	9
2.1 Magnetic Resonance Imaging - Possible Interactions	9
2.1.1 Problematic interactions	9
2.1.2 Fields in the MRI scanner	9
2.1.3 Force of the static and gradient magnetic fields	9
2.1.4 Interactions resulting from induced AC voltages	9
2.1.5 Thermal interactions	10
2.1.6 Image interference and artifacts	10
2.2 Warnings	10
2.2.1 Preliminary notes	10
2.2.2 Warnings	10
3. MRI Conditions for Use	11
3.1 Patient Pre-MRI Conditions	11
3.2 MRI Scanner Limitations	11
3.3 Restrictions during the MR Scan	11
4. MRI Examination	13
4.1 Preliminary Examination	13
4.1.1 Cooperation between specialists	13
4.1.2 Checking the suitability of the patient and the implanted system	13
4.1.3 Performing an MR scan and programming the MRI Mode	13
4.2 MRI Examination	15
4.2.1 Prerequisites	15
4.2.2 Basic conditions and restrictions	15
4.2.3 Patient monitoring during the MR scan	15
4.2.4 Completion of the examination	15

Table of Contents

ProMRI® System Technical Manual

- 5. Post MR Scan Requirements.....17**
 - 5.1 Follow-Up Procedure.....17
- 6. Clinical Study.....19**
 - 6.1 ProMRI Phase B Clinical Study.....19
 - 6.1.1 Primary Objectives19
 - 6.1.2 Methods.....19
 - 6.1.3 Results19
 - 6.2 ProMRI Phase C Clinical Study.....29
 - 6.2.1 Primary Objectives29
 - 6.2.2 Methods.....29
 - 6.2.3 Results29
- 7. Adverse Events.....37**
 - 7.1 Observed Adverse Events.....37
 - 7.1.1 ProMRI Phase B.....37
 - 7.1.2 ProMRI Phase C38

1. Basic Information

1.1 About this Manual

1.1.1 Subject of this manual

This manual provides information about the safe application of an MR scan on patients with a ProMRI® system, which consists of a ProMRI® device and the leads listed in Section 1.4.

This manual will focus on the MRI conditions and safety measures that are to be adhered to before and during an MR scan using the ProMRI® system.

1.1.2 What this manual doesn't include

Correct and safe use of the ProMRI® system components is described in the technical manuals provided with the products and is not a subject of this manual.

Correct and safe use of an MRI scanner is not described in this manual unless directly related to the ProMRI® System.

1.2 Target Group

1.2.1 Cooperation between professionals from two areas of expertise

This manual is intended for physicians and medical staff who have the knowledge and experience required to prepare and perform MR scans on patients with a pacemaker, ICD, or CRT-D.

Preparation and application of an MR scan on a CRM patient requires close cooperation between a cardiology professional, as a specialist for the device system, and a radiology professional, as a specialist for the MR scan.

The following sections describe the tasks that each of these specialists is responsible for.

1.2.2 Knowledge required by the cardiology specialist

A cardiology professional is required to select and/or approve the patient for the MR scan. Additionally, they must test the device system before the exam, program the device to the MRI mode, ensure its functionality after the exam, and program it back to the mode that was active before the MR scan.

The cardiology professional should be knowledgeable in the following areas and subjects:

- Performing pacemaker/ICD/CRT-D therapy
- Handling the BIOTRONIK programmer and especially the following activities:
 - Interrogating the active device
 - Performing follow-up
 - Updating parameters, including programming the device to the MRI Mode
- All associated risks, possible side effects, and the appropriate safety and therapy measures

1.2.3 Knowledge required by the radiology specialist

The radiology professional is responsible for the successful performance of the MR scan for the purposes of the desired diagnosis. Additionally, the radiology professional is also responsible for ensuring that the restrictive conditions, for which the MRI is to be performed, are observed both before and during the MR scan.

The radiology professional should be knowledgeable in the following areas and subjects:

- Handling MRI scanners
- Preparation, performance and analysis of MR scans

1.3 Active Device and Lead

The respective intended use of the ProMRI pacemakers, ICDs, CRT-Ds, and lead(s) applies to use of the device system.

NOTE:

The technical manuals for the components of the ProMRI system are to be observed.

1.3.1 Patient selection, MRI indication

Before a patient with an MR conditional device system is selected for an MR scan, the following issues must be resolved:

- There must be a clear indication for the MR scan.
- Risk/benefit analysis
- All of the exclusion criteria listed in this technical manual have been taken into consideration.
- The described restrictions and conditions for the MR scan are to be observed at all times.

1.3.2 Intended use

If particular MRI conditions are fulfilled, MR scans can now be conducted on patients with a combination of a BIOTRONIK active device and lead that has been tested for this purpose.

1.3.3 Residual risk

The expected risks and hazards are minimized by the measures performed in this manual. Nevertheless, a residual risk remains.

1.4 ProMRI® System

The ProMRI® System consists of the following BIOTRONIK legally marketed devices listed below:

Pacemakers				
Leads	Bradycardia Devices			
	Eluna DR-T (394 969)	Eluna SR-T (394 971)	Entovis DR-T (371 992)	Entovis SR-T (371 994)
	Eluna DR (394 970)	Eluna SR (394 972)	Entovis DR (371 991)	Entovis SR (371 993)
Setrox S 53 (350 974)	ProMRI®	ProMRI®	ProMRI®	ProMRI®
Safio S 53 (370 945)				
Setrox S 60 (350 975)				
Safio S 60 (370 946)				

Continued >

Pacemakers				
Leads	Bradycardia Devices			
	Eluna DR-T (394 969)	Eluna SR-T (394 971)	Entovis DR-T (371 992)	Entovis SR-T (371 994)
	Eluna DR (394 970)	Eluna SR (394 972)	Entovis DR (371 991)	Entovis SR (371 993)
Siello S 45 (362 700)	ProMRI®	ProMRI®	ProMRI®	ProMRI®
Siello S 53 (362 701)				
Siello S 60 (362 702)				
Solia S 45 (377 176)				
Solia S 53 (377 177)				
Solia S 60 (377 179)				

DX ICDs			
Leads	Tachycardia Devices		
	Iperia 7 VR-T DX (DF-1) (393 032)	Inventra 7 VR-T DX (DF-1) (399 436)	Iforia 7 VR-T DX (DF-1) (390 090, 390 092, 390 093, 390 095)
Lincox ^{smart} S DX 65/15 (365 500)	ProMRI®	ProMRI®	ProMRI®
Lincox ^{smart} S DX 65/17 (365 501)			

Chapter 1 Basic Information

DR-T DF-1 ICDs		Tachycardia Devices		
Leads		Iperia 7 DR-T (DF-1) (392 409)	Inventra 7 DR-T (DF-1) (399 430)	Iforia 7 DR-T (DF-1) (390 064, 390 066, 390 067, 390 069)
Atrial	Setrox S 53 (350 974)	ProMRI®	ProMRI®	ProMRI®
	Safio S 53 (370 945)			
	Siello S 45 (362 700)			
	Siello S 53 (362 701)			
	Siello S 60 (362 702)			
	Solia S 45 (377 176)			
	Solia S 53 (377 177)			
	Solia S 60 (377 179)			
ICD	Lincox ^{smart} S 65 (369 818)	ProMRI®	ProMRI®	ProMRI®
	Lincox ^{smart} SD 65/18 (359 067)			

DR-T DF4 ICDs		Tachycardia Devices		
Leads		Iperia 7 DR-T (DF4) (392 423)	Inventra 7 DR-T (DF4) (399 428)	Iforia 7 DR-T (DF4) (390 072, 390 075)
Atrial	Setrox S 53 (350 974)	ProMRI®	ProMRI®	ProMRI®
	Safio S 53 (370 945)			
	Siello S 45 (362 700)			
	Siello S 53 (362 701)			
	Siello S 60 (362 702)			
	Solia S 45 (377 176)			
	Solia S 53 (377 177)			
	Solia S 60 (377 179)			
ICD	Protego S 65 (379 969, 394 099)	ProMRI®	ProMRI®	ProMRI®
	Protego S 75 (379 968, 394 100)			
	Protego SD 65/16 (399 409, 399 414)			
	Protego SD 65/18 (399 410, 399 415)			
	Protego SD 75/18 (399 411, 399 416)			

Chapter 1 Basic Information

HF-T DF-1 ICDs			
Leads		Tachycardia Devices	
		Iperia 7 HF-T (DF-1) (393 007, 399 426, 393 008, 399 427)	Inventra 7 HF-T (DF-1) (393 019, 399 423)
Atrial	Setrox S 53 (350 974)	ProMRI®	ProMRI®
	Safio S 53 (370 945)		
	Siello S 45 (362 700)		
	Siello S 53 (362 701)		
	Siello S 60 (362 702)		
	Solia S 45 (377 176)		
	Solia S 53 (377 177)		
	Solia S 60 (377 179)		
ICD	Linov ^{smart} S 65 (369 818)	ProMRI®	ProMRI®
	Linov ^{smart} SD 65/18 (359 067)		
LV	Corox OTW 75 BP (381 487, 354 805)	ProMRI®	ProMRI®
	Corox OTW 85 BP (381 488, 354 807)		
	Corox OTW-L 75 BP (381 492, 368 345)		
	Corox OTW-L 85 BP (381 491, 368 346)		
	Corox OTW-S 75 BP (381 489, 355 148)		
	Corox OTW-S 85 BP (381 490, 355 149)		

HF-T DF4 ICDs			
Leads		Tachycardia Devices	
		Iperia 7 HF-T (DF4) (393 009, 399 424, 393 010, 399 425)	Inventra 7 HF-T (DF4) (393 020, 399 422)
Atrial	Setrox S 53 (350 974)		
	Safio S 53 (370 945)		
	Siello S 45 (362 700)		
	Siello S 53 (362 701)		
	Siello S 60 (362 702)		
	Solia S 45 (377 176)		
	Solia S 53 (377 177)		
	Solia S 60 (377 179)		
ICD	Protego S 65 (379 969, 394 099)	ProMRI®	ProMRI®
	Protego S 75 (379 968, 394 100)		
	Protego SD 65/16 (399 409, 399 414)		
	Protego SD 65/18 (399 410, 399 415)		
	Protego SD 75/18 (399 411, 399 416)		
LV	Corox OTW 75 BP (381 487, 354 805)		
	Corox OTW 85 BP (381 488, 354 807)		
	Corox OTW-L 75 BP (381 492, 368 345)		
	Corox OTW-L 85 BP (381 491, 368 346)		
	Corox OTW-S 75 BP (381 489, 355 148)		
	Corox OTW-S 85 BP (381 490, 355 149)		

The included BIOTRONIK BS DF-1 blind plugs are approved as MR conditional. BIOTRONIK's BS IS-1 blind plugs are certified as MR conditional when used in the LV connector port of the respective device.

Since the ProMRI® pacemakers, ICDs, CRT-Ds, and the leads are sold independently of each other, this manual informs the user about the MRI conditions for use that are to be observed.

2. Safety Warnings

2.1 Magnetic Resonance Imaging - Possible Interactions

2.1.1 Problematic interactions

Significant mechanisms which can lead to problematic interactions with device systems are described here. Therefore, MR scans are generally contraindicated for cardiac pacemaker/ICD/CRT-D patients. BIOTRONIK has developed the ProMRI® System, which minimizes the effects listed below on the device system and the patients.

2.1.2 Fields in the MRI scanner

The following three types of fields are generated in an MR scan:

- **Static magnetic field**
This is a consistently strong, rectified magnetic field, which is constantly emitted in the MRI scanner and its immediate surroundings, even if no scan is being performed.
- **Gradient magnetic fields**
These are low-frequency pulsed magnetic fields with a relatively low amplitude.
During the MR scan, the patient is exposed to three vertical gradient magnetic fields that are facing towards each other.
- **HF field (high frequency field)**
This is a high frequency electromagnetic field which activates the protons on their resonance frequency. It is switched on several times during the imaging process but only for very short periods.
The HF field is created by so-called emitting coils, which also serve as reception coils.
Differentiation is made between the emitting coils (body coils) integrated in the MRI scanner with the addition of optional local emitting coils (e.g. head coil with transmitting function).

2.1.3 Force of the static and gradient magnetic fields

Implanted ferromagnetic materials are subject to the force of these magnetic fields. This means that implanted devices can subject the surrounding tissue to pressure, tensile force or vibrations. The construction and choice of material in the MR conditional devices and compliance with the specified conditions serve to reduce these stresses to an acceptable minimum.

2.1.4 Interactions resulting from induced AC voltages

Gradient magnetic fields and electromagnetic high frequency fields can induce electrical AC voltages in metallic devices. In some cases, these electrical energies can result in undesirable pacing or have a negative impact on the implanted device.

Constructive measures on the MR conditional devices and the restrictive prerequisites for the arrangement and conduction of the MR scan reduce the probability of occurrence and strength of this effect. However, this effect cannot be entirely excluded.

Among other things, corresponding emergency precautions in this case have to be taken.

2.1.5 Thermal interactions

High-frequency electromagnetic fields induce electric voltages in the lead, which cause current conduction through the lead and the tissue electrically connected to the lead. This flow of current in turn causes warming at the electrical points of contact between the lead and the tissue, which can result in thermal damage to the surrounding tissue. This thermal tissue damage can be temporary or lasting and can cause deterioration of the lead's pacing and sensing functions.

Gradient magnetic fields can cause warming of the device housing, which can lead to thermal exposure and damage to the surrounding tissue.

Due to the constructive composition of the MR conditional devices and the compliance with the tested conditions and restrictions for the MR scan, these thermal effects are kept to a tolerable measure.

2.1.6 Image interference and artifacts

Not only can the MR scan have undesirable effects on the patient or the device system, but the implanted devices can also have a negative impact on the MR scan.

If the devices are outside the scanning area, they can cause slight image distortion and interference.

If a device is within the area shown by the MRI scanner, then artifacts, distortion and interference are probable. Consider this when selecting the image calculation parameters and the depicted area.

2.2 Warnings

2.2.1 Preliminary notes

Please refer to the technical manuals for the ProMRI pacemakers, ICDs, CRT-Ds, and leads. This manual only deals with aspects that are relevant within the MR scan context.

This manual does not deal with the contraindications of MRI applications, which do not result from interactions with a device system.

2.2.2 Warnings

- An MR scan on a device system patient is always contraindicated for device systems that have not been identified as MR conditional by BIOTRONIK and have not been approved for MRI applications by the FDA.
- An MR scan on a patient with an MR conditional device system is also contraindicated when any of the conditions listed in the MRI Conditions for Use ([Section 3](#)) are not adhered to.
- The MR Conditional device system must be programmed to the MRI Mode prior to the MR scan.

3. MRI Conditions for Use

3.1 Patient Pre-MRI Conditions

The following requirements must always be fulfilled in order to perform an MR scan using BIOTRONIK's ProMRI® System:

- The device system consists of a pacemaker, ICD, or CRT-D with the respective leads that are separately labeled MR conditional and, when combined, constitutes an MR conditional device system (See [Section 1.4](#)).
- There are no other active or abandoned cardiac implants (e.g., lead extensions, lead adapters or abandoned leads) in the patient's body.
- Other active or passive implants are permitted if they are identified as MR conditional by the manufacturer.

Note: An MRI scan is permitted only if the product-specific conditions are met for all implants and if no metal implantable device longer than 5 cm is in the vicinity of a BIOTRONIK lead within a distance of less than 4 cm.

- The leads have been implanted for at least 6 weeks.
- The device system is implanted pectorally.
- The measured pacing threshold is not above 2.0 V at 0.4 ms pulse width.
- The device system should be functioning normally prior to an MRI.
- The battery status is neither ERI nor EOS.
- The device is programmed to an MRI mode immediately before the MR scan.

3.2 MRI Scanner Limitations

The MRI scanner has to meet the following conditions:

- Use of a clinical MRI system with a cylindrical bore and a static magnetic field strength of 1.5 Tesla.
- The slew rate of the MRI scanner's gradient fields should not exceed 200 T/m/s per axis.

3.3 Restrictions during the MR Scan

The following conditions must be met during the MR scan:

- The MR scan should be performed with the patient in supine position.
- The mean specific absorption rate (SAR) for the whole body displayed by the MR scanner must not exceed 2.0 W/kg.
- The head absorption rate displayed by the MRI scanner must not exceed 3.2 W/kg.
- Emergency equipment for resuscitation must be kept at hand and properly certified staff must be available.
- Continuously monitor the patient's condition during the entire MR scan using at least one of the following parameters: blood oxygen saturation, blood pressure or ECG.

Note: The ECG function integrated in the MRI scanner is often not permitted for patient monitoring. Therefore, only use devices which are permitted for patient monitoring in an MRI environment.

4. MRI Examination

4.1 Preliminary Examination

4.1.1 Cooperation between specialists

Preparation and conduction of an MR scan on a patient with the ProMRI® System requires close cooperation between a specialist for the device system and a specialist for the MRI technology and MR scan.

One of these specialists has to perform the steps described in the following for preparation of the MR scan, the patient and his or her device system.

The person responsible for each task depends on the activity or context of the scan.

4.1.2 Checking the suitability of the patient and the implanted system

Cardiology and radiology professionals are required for this step, proceed as follows:

Step	Action
1	Check and ensure that all requirements pertaining to the patient and the device system described in the ProMRI® System (Section 1.4) and MRI Conditions for Use (Section 3) are met.
2	Make sure the technical and clinical basic conditions for the MR scan can be met and that the necessary preparations have been made.

4.1.3 Performing an MR scan and programming the MRI Mode

Once the conditions for an MR scan have been clarified, preliminary examination and programming to an MRI mode by the cardiologist are the final and definitive preparation measures; proceed as follows:

Step	Action
1	Apply the programming head of the programmer to the chest and interrogate the device.
2	Perform full follow-up and check the following preconditions for the MR scan: <ul style="list-style-type: none"> • Normal device functionality • Battery status is neither ERI nor EOS • Pacing threshold: max. 2.0 V / 0.4 ms
Programming the MRI Mode	
3	Open the MRI program using one of the following options: <ul style="list-style-type: none"> • Select Follow-up → MRI. • Select Parameters → Bradycardia → Program sets → Show MRI program • Select Parameters → Bradycardia → MRI Mode (for Eluna and ICDs only)

Step	Action
4	<p>Read the preliminary and basic conditions in the MRI checklist window precisely and activate the checkbox I accept the conditions for MRI examinations.</p> <div data-bbox="446 336 1234 919" style="border: 1px solid black; padding: 10px; margin: 10px auto; width: fit-content;"> <p>MRI checklist</p> <p>Check device and leads</p> <ul style="list-style-type: none"> - ProMRI leads have been implanted for at least 6 weeks. - ProMRI device has been implanted in the pectoral region. - No other active or abandoned cardiac devices are present. - Other MR conditional implants may be allowed as specified in the ProMRI user manual. - Pacing threshold (< 2.0 V/0.4 ms) and lead impedance ranges as specified in the ProMRI user manual. <p>Radiological considerations</p> <ul style="list-style-type: none"> - Standard cylindrical scanner architecture required - Check ProMRI manual for maximum field strength allowed - Continuous patient monitoring required during MR scan - Observe specific conditions for MR Conditional devices (SAR, scan zone, field strength, slew rate ...) - After MRI scan, restore previously programmed parameters, program and confirm settings. <p>MRI program ON</p> <p>MRI mode DOO</p> <p>Basic rate [bpm] 90</p> <p>Ventricular pacing RV</p> <p><input checked="" type="checkbox"/> I accept the conditions for MRI examinations</p> <ul style="list-style-type: none"> - Please select an MRI mode. After clicking 'OK' the MRI parameters will be displayed. - Selecting 'Program' on the next screen will program the device accordingly. - Any parameter change will result in the loss of the MRI program. <p style="text-align: center;"> Print Help OK Cancel </p> </div> <p>The software will not permit further programming of an MRI mode without this confirmation.</p>

WARNING

Health risk to patients due to limited device function.

Continuous cardiac monitoring of the patient must be ensured until the device system's full functionality is restored in the follow-up examination.

5	<p>Activate one of the possible MRI modes.</p> <ul style="list-style-type: none"> • The OFF setting is recommended for patients not dependent on their pacemaker or ICD. • An asynchronous mode (DOO, AOO, VOO) is recommended for pacing-dependent patients depending on the particular indication. • DOO-BiV or VOO-BiV (available in CRT-Ds only) is recommended for pacing-dependent patients with a triple-chamber device. <p>Select the basic rate (for Eluna, ICDs, and CRT-Ds only). (For Entovis, the basic rate is set to 80 bpm.)</p> <p>The following parameters are set for the implant:</p> <ul style="list-style-type: none"> • Pulse amplitude: 4.8 V (pacemakers); 5.0 V (ICDs/CRT-Ds) • Pulse width: 1.0 ms • All automatic functions are deactivated. • Home Monitoring remains as programmed (For Entovis only, Home Monitoring is deactivated.) • The magnet response is set to SYNC (synchronous). • ICD therapy is inactive. • For CRT-Ds, programmed settings for the LV lead are maintained if certain bipolar restrictions are met (i.e., LV pacing polarity must be programmed to LV tip → LV ring or LV ring → LV tip).
---	---

NOTE:

You will leave the MRI mode if you change parameters prior to programming the device.

Step	Action
6	Transmit the MRI mode to the implant.
<p>NOTE:</p> <p>When programming the MRI mode in Entovis or Iforia, the original device settings are saved in the programmer. These settings can be accessed again during the follow-up examination after completion of the MR scan, which simplifies restoration of the status from before the MR scan. The same programmer must be used as for the preliminary examination.</p> <p>For Eluna pacemakers, Iperia/Inventra ICDs, or the CRT-Ds, the original settings are stored in the device when the MRI mode is programmed. Therefore, the permanent program can be restored even when the device is interrogated by a different programmer after the MR scan.</p>	
7	Print and document follow-up data (print report).
8	<p>Finish the preliminary examination of the patient.</p> <p>Several device functions may be deactivated in the MRI mode. Therefore, make sure that the patient can be scheduled for a follow-up immediately after completion of the MR scan in order to reprogram the device back to the permanent parameters as defined by the patient's physician.</p>

4.2 MRI Examination

4.2.1 Prerequisites

The following conditions have to be met:

- The contraindications listed in the respective sections as well as the required MRI Conditions for Use are taken into consideration.
- The patient is previously examined by a cardiology professional.
- The device is programmed to the MRI Mode, which is suitable for an MR scan.
- The technical and organizational conditions are met to be able to observe the restrictions and safety measures required during the MR scan.
- Emergency equipment for resuscitation (including specialist staff certified to use it) is available.

4.2.2 Basic conditions and restrictions

The MRI Conditions for Use ([Section 3](#)) have to be met during an MR scan and the device must be programmed to MRI Mode prior to the MR scan.

4.2.3 Patient monitoring during the MR scan

The patient should be continuously monitored during the entire MR scan, including maintaining visual and verbal contact with the patient and monitoring of blood oxygen saturation, blood pressure or ECG. Emergency equipment for resuscitation must be kept at hand and properly certified staff must be available.

If the patient exhibits signs of discomfort (i.e., warming is noted) or hemodynamic function appears to be compromised at any point during the scan, discontinue the scan and remove the patient from the MRI scanner.

4.2.4 Completion of the examination

After completing the MR scan, make sure the patient is monitored by a cardiology professional who performs the required follow-up cardiology examination and reprogramming of the device system.

5. Post MR Scan Requirements

After the MR scan, the patient should undergo a follow-up device interrogation. This is necessary for the patient's safety for two reasons:

- To reprogram the device back into the original parameters.
- To assess the device system for potential adverse effects caused by the MR scan.

NOTE:

The device parameters during activation of the MRI mode are maintained until the MRI program is set to Off after the MR scan.

5.1 Follow-Up Procedure

After an MR scan, the follow-up procedure should be performed:

1. Apply the programming head.
2. Interrogate the device.
3. Reactivate the program that was effective prior to programming the MRI mode.
4. Send the reactivated program to the device.
5. Perform a complete follow-up.
6. If necessary, perform further examinations.
7. Print and document follow-up data (print report).
8. Finish the follow-up for the patient.

6. Clinical Study

6.1 ProMRI Phase B Clinical Study

6.1.1 Primary Objectives

This clinical investigation was designed to demonstrate the clinical safety of the ProMRI Pacemaker System when used under specific magnetic resonance imaging (MRI) conditions for full body MRI scan. The investigation included five primary endpoints, which condense into three main objectives:

- Primary Endpoint 1 – Evaluation of serious adverse device effect (SADE) rate related to the implanted pacing system and MRI procedure
- Primary Endpoints 2 & 3 – Evaluation of atrial and ventricular lead pacing threshold increases
- Primary Endpoints 4 & 5 – Evaluation of P-wave and R-wave sensing attenuation

6.1.2 Methods

The study enrolled subjects implanted with an Entovis family pacemaker (SR-T, DR-T) and one or two Setrox S 53 or 60 leads, and were willing to undergo an MRI scan.

The patients selected for participation were from the investigator's general patient population meeting the indications for use of the Entovis family pacemaker system. To qualify for enrollment, subjects were required to have measurable pacing thresholds ≤ 2.0 V @ 0.4 ms and could not be implanted with other non-MRI compatible devices. Patients received a baseline evaluation at least seven days prior to the MRI procedure, at which time the pacemaker was tested and programmed to an MRI mode before the MRI, then tested and reprogrammed to the original pacing mode post-MRI. The study required a cardiac or thoracic spine MRI scan.

Patients were enrolled post-implant, underwent an MRI procedure and testing, and were followed at one and three months post-MRI. During follow-up visits, a device interrogation was completed and the investigator determined if the MRI scan had any long-term effects on the function of the pacemaker system.

6.1.3 Results

A total of 244 subjects were provisionally enrolled and 216 subjects were fully enrolled at 31 sites as of August 12, 2014. The cumulative implant duration of the 216 fully enrolled subjects at baseline and MRI procedure was 96.5 years (average implant duration of 0.45 ± 0.32 years) and 111.6 years (average implant duration of 0.52 ± 0.33), respectively. The patient follow-up compliance rate was 98.8% out of 341 required follow-ups. Endpoint data is provided for the Per Protocol (PP) and Intention-to-treat (ITT) Populations. The ITT population includes all subjects programmed to MRI mode with endpoint data from follow-up or Home Monitoring. At the time of data analysis, 199 had completed their 1-month follow-up. An additional four subjects had a missed 1-month follow-up, but are included in the ITT endpoint analysis using their Home Monitoring data. The average subject is a 68 year old male who weighs 191 pounds and is 68 inches in height.

Primary Endpoint 1

The purpose of Primary Endpoint 1 was to evaluate the rate of Serious Adverse Device Effects related or possibly related to the implanted pacing system and the MRI procedure. Only SADEs that were pacing system and MRI related or possibly related, as adjudicated by the independent Data Monitoring Committee, were taken into account for calculation of the SADE rate.

Analysis

The Data Monitoring Committee (DMC) adjudicated 25 events reported by the investigators. Of the 25 events reported by the investigators, one was adjudicated as related or possibly related to both the implanted pacing system and the MRI procedure resulting in an SADE-free rate of 99.5% (202/203), $p < 0.001$, 95% CI: (97.3%, 100.0%). Two additional subjects experienced events that were adjudicated as possibly related to the MRI procedure only. Accounting for the two additional events that were possibly related to the MRI, the estimated freedom from Adverse Events was 98.5% (200/203), $p < 0.001$, 95% CI: (95.7%, 99.7%). A rejection of the null hypothesis indicates that the SADE-free rate possibly related to the implanted pacing system and the MRI procedure is greater than 90% at 1-month post-MRI and the endpoint is met.

Primary Endpoints 2 & 3

The purpose of Primary Endpoints 2 and 3 was to evaluate the percentage of atrial and ventricular pacing leads with a pacing threshold increase between the pre-MRI and 1-month post-MRI follow-up. The threshold behavior of the lead is defined as a success if the increase is not larger than 0.5 V. Table 1, Figure 1, and Figure 2 display the differences in atrial pacing thresholds and Table 2, Figure 3, and Figure 4 display the difference in ventricular pacing thresholds.

Table 1: Atrial Pacing Threshold

	Results	P-value*
Intention-to-treat: Includes Values Imputed from Home Monitoring (N=189)		
Difference in Atrial Pacing Threshold (V)		
Mean ± SD (N)	0.00 ± 0.13	
Minimum, Median, Maximum	-0.70, 0.00, 0.40	
Leads Meeting Success Criteria (%)	189 (100.0%)	P < 0.001
95% Confidence Interval	(98.1%, 100.0%)	
Per Protocol (N=178)		
Difference in Atrial Pacing Threshold (V)		
Mean ± SD (N)	0.00 ± 0.12	
Minimum, Median, Maximum	-0.70, 0.00, 0.40	
Leads Meeting Success Criteria (%)	178 (100.0%)	P < 0.001
95% Confidence Interval	(97.9%, 100.0%)	

*Exact binomial test (1-sided) for comparison to 95%

Figure 1: Histogram of PPP Atrial Pacing Threshold Differences (One-Month – Pre-MRI)





Figure 2: Histogram of ITT Atrial Pacing Threshold Differences (One-Month – Pre-MRI)



Atrial Analysis

The mean threshold increase for the PP and ITT populations was 0.00 ± 0.12 V and 0.00 ± 0.13 V, respectively. Of 178 total subjects in the PP population and the 189 total subjects with data in the ITT population, all had a change in atrial pacing threshold of less than or equal to 0.5V between one-month post-MRI and pre-MRI. A rejection of the null hypothesis (p-values: PP – <0.001 , ITT – <0.001) indicates that the proportion of atrial pacing threshold success is greater than 95% and Primary Endpoint 2 is met.

Table 2: Ventricular Pacing Threshold

	Results	P-value*
Intention-to-treat: Includes Values Imputed from Home Monitoring (N=199)		
Difference in Ventricular Pacing Threshold (V) Mean ± SD (N) Minimum, Median, Maximum	0.00 ± 0.10 -0.30, 0.00, 0.30	
Leads Meeting Success Criteria (%) 95% Confidence Interval	199 (100.0%) (98.2%, 100.0%)	P < 0.001
Per Protocol (N=189)		
Difference in Ventricular Pacing Threshold (V) Mean ± SD (N) Minimum, Median, Maximum	0.01 ± 0.10 -0.20, 0.00, 0.30	
Leads Meeting Success Criteria (%) 95% Confidence Interval	189 (100.0%) (98.1%, 100.0%)	P < 0.001

*Exact binomial test (1-sided) for comparison to 95%

Figure 3: Histogram of PPP Ventricular Pacing Threshold Differences (One-Month – Pre-MRI)





Figure 4: Histogram of ITT Ventricular Pacing Threshold Differences (One-Month – Pre-MRI)



Ventricular Analysis

The mean threshold increase for the PP populations was 0.01 ± 0.10 and for the ITT population was 0.00 ± 0.10 . Of 189 total subjects in the PP population and the 199 total subjects with data in the ITT population, all had a change in ventricular pacing threshold of less than or equal to 0.5V between one-month post-MRI and pre-MRI, respectively. A rejection of the null hypothesis (p -values: PP – <0.001 and ITT – <0.001) indicates that the proportion of ventricular pacing threshold success is greater than 95% and Primary Endpoint 3 is met.

Primary Endpoints 4 & 5

The purpose of Primary Endpoints 4 and 5 was to evaluate the percentage of subjects who experienced P-wave and R-wave attenuation between the pre-MRI and 1-month post-MRI follow-up. Sensing amplitude attenuation was defined as either a P-wave or R-wave amplitude decrease (between pre-MRI and one month follow-up) exceeding 50% or an amplitude at the one month follow-up of less than 1.5 mV and 5.0 mV in the atrium and ventricle, respectively. Table 3 and Table 4 display differences in atrial and ventricular sensing amplitudes over this period of time. Figure 5 and Figure 6 display the PP atrial sensing amplitude ratios with endpoint boundary conditions. Figure 7 and Figure 8 display the ITT sensing amplitude ratios with endpoint boundary conditions.

Table 3: P-Wave Sensing Attenuation

	Results	P-value*
Intention-to-treat): Includes Values Imputed from Home Monitoring (N=189)		
Difference in P-wave Amplitude (mV) Mean ± SD (N) Minimum, Median, Maximum	-0.14 ± 0.88 -2.7, -0.1, 2.8	
Leads Meeting Success Criteria (%) 95% Confidence Interval	186 (98.4%) (95.4%, 99.7%)	P < 0.001
Per Protocol (N=157)		
Difference in P-wave Amplitude (mV) Mean ± SD (N) Minimum, Median, Maximum	-0.21 ± 0.89 -2.7, 0.0, 2.8	
Leads Meeting Success Criteria (%) 95% Confidence Interval	154 (98.1%) (94.5%, 99.6%)	P < 0.001

*Exact binomial test (1-sided) for comparison to 90%

Figure 5: Atrial PPP Sensing Amplitude Ratio with Boundary Conditions





Figure 6: Atrial ITT Sensing Amplitude Ratio with Boundary Conditions



Atrial Analysis

Of 157 total subjects in the PP population and 189 total subjects with data in the ITT population, 154 (98.1%, $p < 0.001$) and 186 (98.4%, $p < 0.001$) met the endpoint for attenuation-free P-wave sensing, respectively. A rejection of the null hypothesis indicates that the P-wave attenuation free rate is greater than 90% and Primary Endpoint 4 is met.

Table 4: R-Wave Sensing Attenuation

	Results	P-value*
Intention-to-treat: Includes Values Imputed from Home Monitoring (N=199)		
Difference in R-wave Amplitude (mV)		
Mean \pm SD (N)	-0.24 \pm 1.39	
Minimum, Median, Maximum	-6.1, -0.1, 4.5	
Leads Meeting Success Criteria (%)	199 (100.0%)	P < 0.001
95% Confidence Interval	(98.2%, 100.0%)	
Per Protocol (N=173)		
Difference in R-wave Amplitude (mV)		
Mean \pm SD (N)	-0.26 \pm 1.37	
Minimum, Median, Maximum	-6.1, -0.1, 4.5	
Leads Meeting Success Criteria (%)	173 (100.0%)	P < 0.001
95% Confidence Interval	(97.9%, 100.0%)	

*Exact binomial test (1-sided) for comparison to 90%

Figure 7: Ventricular PPP Sensing Amplitude Ratio with Boundary Conditions





Figure 8: Ventricular ITT Sensing Amplitude Ratio with Boundary Conditions



Ventricular Analysis

Of the 173 subjects in the PP population and 199 subjects with data in the ITT population, all met the endpoint for attenuation-free R-wave sensing. A rejection of the null hypothesis (p-values: PP – <0.001 and ITT – <0.001) indicates that the R-wave attenuation free rate is greater than 90% and Primary Endpoint 5 was met.

Additional Data of Interest: Multiple MRI Scans

Subjects both with and without indications for MRI scans could be included in the study. At least one MRI scan per subject was required in the course of the study.

Table 5 provides the number of Phase B subjects who have undergone multiple MRI scans.

Table 5: Subjects with Multiple MRI Scans

MRI Scan Type	Number of Subjects
Subjects with 3 MRI Scans: <ul style="list-style-type: none"> Protocol required Phase A MRI Scans (Head and Pelvic MRI Scan) Protocol required Phase B MRI Scan (Thoracic Spine or Cardiac MRI Scan) 	76
Subjects with 4 MRI Scans: <ul style="list-style-type: none"> Protocol required Phase A MRI Scans (Head and Pelvic MRI Scan) Protocol required Phase B MRI Scan (Thoracic Spine or Cardiac MRI Scan) Clinically Indicated MRI Scan (Phase A or Phase B) 	5

There were 81 subjects of the 216 total subjects who completed the Phase B protocol that contributed data for the following analyses.

Table 6 provides a summary of the atrial and ventricular pacing threshold difference, P-wave and R-wave sensing amplitude difference and atrial and ventricular impedance difference between the post-MRI Phase B procedure and pre-MRI Phase B procedure, as well as between the post-MRI Phase B procedure and pre-MRI Phase A procedure.

Table 6: Atrial Pacing Threshold Differences

N Mean ± SD Range	Difference between post-MRI Phase B and pre-MRI Phase B procedures	Difference between post-MRI Phase B and pre-MRI Phase A procedures
Atrial Pacing Threshold (V)	74 -0.01 +/- 0.10 -0.3 to 0.3	76 -0.01 +/- 0.17 -0.3 to 0.6
Ventricular Pacing Threshold (V)	78 0.02 +/- 0.06 -0.1 to 0.1	78 0.03 +/- 0.15 -0.4 to 0.4
P-Wave Sensing Amplitude (mV)	76 -0.14 +/- 0.95 -5.0 to 2.9	76 -0.29 +/- 0.92 -3.9 to 2.3
R-Wave Sensing Amplitude (mV)	79 -0.20 +/- 1.16 -3.0 to 5.6	78 -0.11 +/- 2.03 -4.7 to 5.9
Atrial Pacing Impedance (Ω)	74 -1.0 +/- 13.2 -39 to 39	74 -2.8 +/- 35.0 -98 to 78
Ventricular Pacing Impedance (Ω)	76 -1.5 +/- 10.8 -39 to 20	75 -18.7 +/- 44.2 -253 to 78

Multiple MRI Analysis

None of the subjects with multiple MRIs experienced an endpoint related increase in pacing threshold. Two subjects experienced an endpoint related P-wave sensing attenuation following their Phase B MRI procedure. Both subjects had normal, stable P-wave sensing amplitudes throughout the study as evidenced by Home Monitoring.

Overall Results Summary

The ProMRI® Study is designed to demonstrate the clinical safety of the ProMRI® Pacemaker System when used under specific MRI conditions. A total of 244 subjects were provisionally enrolled at 31 sites as of August 12, 2014.

- The Data Monitoring Committee adjudicated 25 events reported by investigators as hospitalizations or possibly related to the MRI. One adverse event was adjudicated as possibly related to both the implanted pacing system and the MRI procedure, resulting in a SADE-free rate of 99.5% (202/203), $p < 0.001$, 95% CI: (97.3%, 100.0%). A rejection of the null hypothesis indicates that the SADE-free rate possibly related to the implanted pacing system and the MRI procedure is less than or equal to 90% at 1-month post-MRI and the endpoint is met.
- Of 178 per protocol (PP) and 189 intention-to-treat (ITT) subjects, 178 (100.0%) and 189 (100.0%) had a change in atrial pacing threshold of less than or equal to 0.5V between one-month post-MRI and pre-MRI, respectively, resulting in a PP p-value of < 0.001 , 95% CI: (97.9%, 100.0%) and an ITT p-value of < 0.001 , 95% CI: (98.1%, 100.0%). A rejection of the null hypothesis indicates that the proportion of atrial pacing threshold success is greater than 95% and the endpoint is met.
- Of 189 per protocol (PP) and 199 intention-to-treat (ITT) subjects, all had a change in ventricular pacing threshold of less than or equal to 0.5V between one-month post-MRI and pre-MRI resulting in a PP p-value of < 0.001 , 95% CI: (98.1%, 100.0%) and an ITT p-value of < 0.001 , 95% CI: (98.1%, 100.0%). A rejection of the null hypothesis indicates that the proportion of ventricular pacing threshold success is greater than 95% and the endpoint is met.
- Of 157 per protocol (PP) and 189 intention-to-treat (ITT) subjects, 154 (98.1%) and 186 (98.4%) met the endpoint for attenuation-free P-wave sensing, respectively, resulting in a PP p-value of < 0.001 , 95% CI: (94.5%, 99.6%) and an ITT p-value of < 0.001 , 95% CI: (95.4%, 99.7%). A rejection of the null hypothesis indicates that the P-wave attenuation free rate is greater than 90% and the endpoint is met.
- Of 173 per protocol (PP) and 199 intention-to-treat (ITT) subjects, 173 (100.0%) and 199 (100.0%) met the endpoint for attenuation-free R-wave sensing, respectively, resulting in a PP p-value of < 0.001 , 95% CI: (97.9%, 100.0%) and an ITT p-value of < 0.001 , 95% CI: (98.2%, 100.0%). A rejection of the null hypothesis indicates that the R-wave attenuation free rate is greater than 90% and the endpoint is met.

All five primary endpoints were met. The data received and analyzed demonstrates and supports the clinical safety and efficacy of the ProMRI® Pacemaker System when used under specific MRI conditions and without an exclusion zone.

6.2 ProMRI Phase C Clinical Study

6.2.1 Primary Objectives

This clinical investigation was designed to demonstrate the clinical safety of the ProMRI ICD System when used under specific magnetic resonance imaging (MRI) conditions for full body MRI scan. The investigation included three primary endpoints, which condense into three main objectives:

- Primary Endpoint 1 – Evaluation of serious adverse device effect (SADE) rate related to the implanted pacing system and MRI procedure
- Primary Endpoint 2 – Evaluation of ventricular lead pacing threshold increases
- Primary Endpoint 3 – Evaluation of R-wave sensing attenuation

6.2.2 Methods

The study enrolled subjects implanted with an ICD System consisting of an Iforia DR-T and Linox^{smart} S 65 or Linox^{smart} SD 65/18 ICD lead with a Setrox S 53 atrial lead, or Iforia VR-T DX and Linox^{smart} S DX 65/15 or Linox^{smart} S DX 65/17 ICD lead, and were willing to undergo an MRI scan.

The patients selected for participation were from the investigator's general patient population meeting the indications for use of the Iforia ICD system. To qualify for enrollment, subjects were required to have measurable pacing thresholds ≤ 2.0 V @ 0.4 ms and could not be implanted with other non-MRI compatible devices. Patients received a baseline evaluation at least 7 days prior to the MRI procedure, at which time the ICD system was tested and programmed to an MRI mode before the MRI, then tested and reprogrammed to the original programmed parameters post-MRI. For Phase C of the study, the protocol required cardiac or thoracic spine MRI scans.

Patients were enrolled post-implant, underwent an MRI procedure and testing, and were followed at one and three months post-MRI. During follow-up visits, a device interrogation was completed and the investigator determined if the MRI scan had any long-term effects on the function of the ICD system.

6.2.3 Results

A total of 170 subjects were provisionally enrolled at 39 sites. There were 16 subjects that did not meet the MRI criteria or were exited prior to the MRI procedure. The remaining 154 subjects were fully enrolled and had a cumulative implant duration at the time of enrollment of 34.0 years (average implant duration of 0.22 ± 0.14 years). Endpoint data is provided for the Per Protocol (PP) and Intention-to-treat (ITT) Populations. The ITT population includes all subjects programmed to MRI mode with endpoint data from follow-up or Home Monitoring. 154 patients were programmed into the MRI mode at their MRI visit and 151 had completed their 1 month follow-up (2 subjects had a missed visit, one subject exited). The average subject is a 60 year old male who weighs 200.9 pounds and is 68.3 inches in height. The patient follow-up compliance rate was 99.3% out of 303 required follow-ups.

Primary Endpoint 1

The purpose of Primary Endpoint 1 was to evaluate the rate of Serious Adverse Device Effects related or possibly related to the implanted ICD system and the MRI procedure at 1 month post-MRI. Only SADEs that were ICD system and MRI related or possibly related, as adjudicated by the independent Data Monitoring Committee, were taken into account for calculation of the SADE rate.

Primary Endpoint 1 Analysis

The DMC adjudicated 52 Adverse Events reported by the investigators. Of these, none were adjudicated as related or possibly related to both the implanted pacing system and the MRI procedure resulting in an SADE-free rate of 100% (154/154), $p < 0.001$, 95% CI: (97.6%, 100.0%). A rejection of the null hypothesis indicates that the SADE-free rate possibly related to the implanted pacing system and the MRI procedure is greater than 90% at 1 month post-MRI and Primary Endpoint 1 was met.

Primary Endpoints 2

The purpose of Primary Endpoint 2 was to evaluate the percentage of ventricular pacing leads with a pacing threshold increase between the pre-MRI and 1-month post-MRI follow-up. The threshold behavior of the lead is defined as a success if the increase is not larger than 0.5 V. Table 7 displays the ventricular pacing threshold difference between the one-month post-MRI procedure and pre-MRI procedure. Figure 9 displays a histogram of ventricular pacing threshold differences between one-month and pre MRI for the per protocol population and Figure 10 displays a histogram of the data available for ventricular pacing threshold differences for the intention-to-treat population.

Table 7: Ventricular Pacing Threshold

Group	Results	P-value*
Intention-to-treat (ITT): Includes Values Imputed from Home Monitoring (N=154)		
Difference in Ventricular Pacing Threshold (V)		
Mean \pm SD (N)	-0.01 \pm 0.12 (154)	
Minimum, Median, Maximum	-0.40, 0.00, 0.50	
% of Leads Meeting Success Criteria	100.0% (154/154)	P <0.001
95% Confidence Interval	(97.6%, 100.0%)	
Per Protocol (N=147)		
Difference in Ventricular Pacing Threshold (V)		
Mean \pm SD (N)	-0.01 \pm 0.13 (147)	
Minimum, Median, Maximum	-0.40, 0.00, 0.50	
% of Leads Meeting Success Criteria	100.0% (147/147)	P <0.001
95% Confidence Interval	(97.5%, 100.0%)	

*Exact binomial test (1-sided) for comparison to 95%

Figure 9: Histogram of PP Ventricular Pacing Threshold Differences (One-Month – Pre-MRI)





Figure 10: Histogram of ITT Ventricular Pacing Threshold Differences (One-Month – Pre-MRI)



Primary Endpoint 2 Analysis

The mean threshold increase for the PP and ITT population was -0.01 ± 0.13 V and -0.01 ± 0.12 V, respectively. Of 147 total subjects in the PP population and the 154 total subjects with data in the ITT population, 147 (100.0%) and 154 (100.0%) had a change in ventricular pacing threshold of less than or equal to 0.5V between one-month post-MRI and pre-MRI, respectively. A rejection of the null hypothesis (p-values: PP – <0.001 and ITT – <0.001) indicates that the proportion of ventricular pacing threshold success is greater than 95% and Primary Endpoint 2 was met.

Primary Endpoints 3

The purpose of Primary Endpoint 3 was to evaluate the percentage of subjects who experienced R-wave attenuation between the pre-MRI and 1-month post-MRI follow-up. Sensing amplitude attenuation was defined as an R-wave amplitude decrease (between pre- MRI and one month follow-up) exceeding 50% or an amplitude at the one month follow-up of less than 5.0 mV.

Table 8 displays the R-wave sensing amplitude differences between one-month post-MRI procedure and pre-MRI procedure. Figure 11 displays the ventricular sensing amplitude ratio with boundary conditions for the per-protocol population and Figure 12 displays the data available for the intention-to-treat population.

Table 8: R-Wave Sensing Attenuation

	Results	P-value*
Intention-to-treat (ITT): Includes Values Imputed from Home Monitoring (N=153)		
Difference in R-wave Amplitude (mV)		
Mean ± SD (N)	-0.45 ± 1.98 (153)	
Minimum, Median, Maximum	-7.40, 0.00, 6.60	
% of Leads Meeting Success Criteria	99.3% (152/153)	P <0.001
95% Confidence Interval	(96.4%, 100.0%)	
Per Protocol (N=146)		
Difference in R-wave Amplitude m(V)		
Mean ± SD (N)	-0.41 ± 2.01 (146)	
Minimum, Median, Maximum	-7.40, 0.00, 6.60	
% of Leads Meeting Success Criteria	100% (146/146)	P <0.001
95% Confidence Interval	(97.5%, 100.0%)	

*Exact binomial test (1-sided) for comparison to 90%

Figure 11: Ventricular PP Sensing Amplitude Ratio with Boundary Conditions





Figure 12: Ventricular ITT Sensing Amplitude Ratio with Boundary Conditions



Primary Endpoint 3 Analysis

Of 146 total subjects in the PP population and 153 total subjects with data in the ITT population, 146 (100.0%) and 152 (99.3%) met the endpoint for attenuation-free R-wave sensing, respectively. A rejection of the null hypothesis (p-values: per protocol – <0.001 and intent-to-treat – <0.001) indicates that the R-wave attenuation free rate is greater than 90% and Primary Endpoint 3 was met.

Additional Data of Interest: Ventricular Episode Detection Post-MRI

Information about spontaneous VT/VF episodes was collected to determine irregularities in detection or conversion of VT/VF. Investigators were encouraged to program a VT monitoring zone and conduct defibrillation threshold testing (DFT) following the MRI procedure. Table 9 provides a summary of all VT/VF episodes detected from the date of baseline through the pre-MRI procedure and post-MRI through study exit as reported in the EDC system and transmitted by Home Monitoring.

Table 9: Episodes Detected Pre- and Post-MRI

Detection Zone	Baseline through Pre-MRI		Post-MRI Procedure	
	Number of Subjects with Episodes, N (%)	Number of Episodes	Number of Subjects with Episodes, N (%)	Number of Episodes
VT Mon.	8 (5%)	89	27 (18%)	539
VT	2 (1%)	3	11 (7%)	36
VF	5 (3%)	8	11 (7%)	20
Total	15 (10%)	100	49 (32%)	595

Table 10 provides a summary and subsequent success of each VT and VF ventricular therapy sequence. There were 9 episodes in 5 subjects that were excluded due to detections of AF with RVR in the VT/VF therapy zone. An additional 12 were excluded due to no IEGMs available for success determination.

Table 10: Ventricular Therapy Sequences Success Detail

Therapy Sequence	Subjects	Episodes	Successes	Success Rate
ATP	8	10	10	100%
ATP, Shock	2	10	10	100%
Shock	1	2	2	100%
ATP, Spontaneous Conversion	2	3	3	100%
Spontaneous Conversion	7	10	10	100%
All Therapy Sequences	13	35	35	100%

There were no adverse events reported during the study related to inadequate or delayed ICD detection or cardiac arrhythmias. No reports of over or undersensing were noted.

Ventricular Episode Detection Post-MRI Analysis

All episode detections and therapies delivered for VT/VF episodes post-MRI (35/35, 100%) indicate that there was no delayed VT/VF detection due to the MRI procedure. Results demonstrate that the ProMRI® ICD System provides appropriate detection and effective therapy for the treatment of ventricular arrhythmias post-MRI.

Additional Data of Interest: Multiple MRI Scans

Subjects both with and without indications for MRI scans may be included in the study. At least one MRI scan per subject was required in the course of the study.

There were four subjects in Phase C of ProMRI who received a clinically indicated scan in addition to their study MRI scan. Two subjects completed their clinical scan first, and two subjects first underwent their protocol required MRI scan. A clinically indicated post-MRI procedure and three month follow-up procedure were completed on the same day for two subjects.

For the subjects receiving multiple MRI procedures, Table 11, Table 12, and Table 13, respectively, display each ProMRI Phase C subject's ventricular sensing, threshold, and impedance values at the study and clinical MRI procedures and 3-Month follow-up. In this analysis, the clinical and study MRI scans are labeled as MRI1 and MRI2 based on chronological order.

Table 11: ProMRI R-wave Sensing

R-wave Sensing (mV)	MRI 1		MRI 2		3-Month
	Pre	Post	Pre	Post	
Subject A	11.1	11.6	14	12.4	12.4
Subject B	10.7	8.2	8.3	8.6	8.5
Subject C	24.2	24.2	24.2	24.2	24.2
Subject D	24.2	24.2	23.7	24.1	24.2

Table 12: ProMRI Ventricular Threshold

Ventricular Threshold (V)	MRI 1		MRI 2		3-Month
	Pre	Post	Pre	Post	
Subject A	0.4	0.4	0.5	0.5	0.5
Subject B	0.9	0.9	0.9	0.9	1.0
Subject C	1.2	1.1	1.3	1.3	1.3
Subject D	0.5	0.4	0.5	0.5	0.5

Table 13: ProMRI Ventricular Pacing Impedance

Ventricular Impedance (Ω)	MRI 1		MRI 2		3-Month
	Pre	Post	Pre	Post	
Subject A	591	606	606	591	591
Subject B	434	419	434	419	448
Subject C	506	492	463	463	463
Subject D	609	609	609	594	652

Multiple MRI Analysis

None of the subjects with multiple MRIs experienced an endpoint related increase in pacing threshold. Additionally, these subjects did not experience a device system or MRI related adverse event. These data suggest that multiple MRIs do not adversely impact the safety and performance of the system.

Overall Results Summary

The ProMRI study was designed to demonstrate the clinical safety of the ProMRI ICD System when used under specific magnetic resonance imaging (MRI) conditions. A total of 170 subjects were provisionally enrolled at 39 sites; 154 subjects were fully enrolled and 16 subjects did not meet the MRI criteria or were exited prior to the MRI procedure.

- The Data Monitoring Committee adjudicated 52 events reported as serious, hospitalizations, or possibly related to the MRI as reported by the investigators. Of the 52 events reported by the investigators, none were adjudicated as related or possibly related to both the implanted pacing system and the MRI procedure resulting in an SADE-free rate of 100% (154/154), $p < 0.001$, 95% CI: (97.6%, 100.0%). A rejection of the null hypothesis indicates that the SADE-free rate possibly related to the implanted pacing system and the MRI procedure is less than or equal to 90% at 1 month post-MRI and the endpoint is met.
- Of 147 total subjects in the PP population and the 154 total subjects with data in the ITT population, 147 (100.0%) and 154 (100.0%) had a change in ventricular pacing threshold of less than or equal to 0.5V between one-month post-MRI and pre-MRI, respectively. A rejection of the null hypothesis (p -values: per protocol – <0.001 and intent-to-treat – <0.001) indicates that the proportion of ventricular pacing threshold success is greater than 95% and the endpoint is met.
- Of 146 total subjects in the PP population and 153 total subjects with data in the ITT population, 146 (100.0%) and 152 (99.3%) met the endpoint for attenuation-free R-wave sensing, respectively. A rejection of the null hypothesis (p -values: per protocol – <0.001 and intent-to-treat – <0.001) indicates that the R-wave attenuation free rate is greater than 90% and the endpoint is met.

Chapter 6 Clinical Study

ProMRI® System Technical Manual

- All episode detections and therapies delivered for VT/VF episodes post-MRI (35/35, 100%) indicate that there was no delayed VT/VF detection due to the MRI procedure. Results demonstrate that the ProMRI® ICD System provides appropriate detection and effective therapy for the treatment of ventricular arrhythmias post-MRI.

The data received and analyzed demonstrates and supports the clinical safety and efficacy of the ProMRI® ICD System when used under specific MRI conditions without exclusion zone. All three primary endpoints were met with statistical significance. All other additional data of interest were analyzed and support the safety and efficacy of the ProMRI ICD system.

7. Adverse Events

7.1 Observed Adverse Events

7.1.1 ProMRI Phase B

The ProMRI Phase B clinical study data set included 216 enrolled subjects who were programmed into MRI mode with a cumulative number of Subject-Years since enrollment of 60.77.

Adverse events were classified as serious or non-serious. Serious adverse events were defined as events that resulted in a life-threatening illness or injury, resulted in permanent impairment of body structure or function, required in-patient hospitalization, resulted in medical or surgical intervention to prevent life threatening illness or permanent impairment, or led to fetal complications.

Of the 88 adverse events (AEs) reported, there have been 18 serious adverse events (SAEs) in 18 subjects and 70 non-serious adverse events in 66 subjects. A Data Monitoring Committee adjudicated all SAEs and one adverse event was found to be possibly related to both the pacemaker system and the MRI procedure. Two events were adjudicated as possibly related to the MRI procedure.

Table 14: Summary of Serious Adverse Events

Serious Adverse Event	Subjects with Serious Adverse Event, n	Subjects with Serious Adverse Events, %	Serious Adverse Events, n	Serious Adverse Events per Subject-Year
Angina	3	1.4%	3	0.049
Arrhythmia	1	0.5%	1	0.016
Arterial Stenosis	2	0.9%	2	0.033
Gastrointestinal	1	0.5%	1	0.016
Infection	1	0.5%	1	0.016
Medication Related	1	0.5%	1	0.016
Musculoskeletal	2	0.9%	2	0.033
Other	5	2.3%	5	0.082
Syncope/Pre-Syncope	2	0.9%	2	0.033
Total	18	8.3%	18	0.296

Number of Enrolled Subjects = 216, Number of Subject-Years since Enrollment = 60.77

Table 15: Summary of Non-serious Adverse Events

Adverse Event	Subjects with Adverse Event, n	Subjects with Adverse Events, %	Adverse Events, n	Adverse Events per Subject-Year
Arrhythmia	2	0.9%	2	0.033
CAD	1	0.5%	1	0.016
Cancer	1	0.5%	1	0.016
Dermatological	1	0.5%	1	0.016
Fatigue	2	0.9%	2	0.033
Gastrointestinal	2	0.9%	2	0.033
Influenza	1	0.5%	1	0.016
MRI Incidental Finding	30	13.9%	33	0.543
MRI Procedure Discomfort	7	3.2%	7	0.115
Medication Related	4	1.9%	4	0.066
Musculoskeletal	4	1.9%	4	0.066
Other	6	2.8%	7	0.115
Pocket Pain	1	0.5%	1	0.016
Respiratory	1	0.5%	1	0.016
Syncope/Pre-Syncope	2	0.9%	2	0.033
Thrombus	1	0.5%	1	0.016
Total	66	30.6%	70	1.152

Number of Enrolled Subjects = 216, Number of Subject-Years since Enrollment = 60.77

7.1.2 ProMRI Phase C

The ProMRI Phase C clinical study data set included 154 enrolled subjects who were programmed into MRI mode with a cumulative number of Subject-Years since enrollment of 52.3.

Adverse events were classified as serious or non-serious. Serious adverse events were defined as events that resulted in a life-threatening illness or injury, resulted in permanent impairment of body structure or function, required in-patient hospitalization, resulted in medical or surgical intervention to prevent life threatening illness or permanent impairment, or led to fetal complications. All adverse events reported are summarized below as reported by the investigator.

Of the 124 adverse events (AEs) reported, there have been 40 serious adverse events (SAEs) in 23 subjects and 84 non-serious adverse events in 57 subjects. A Data Monitoring Committee adjudicated all SAEs and none of the SAEs were adjudicated as related or possibly related to both the implanted pacing system and the MRI procedure.

Table 16: Summary of Serious Adverse Events

Serious Adverse Event	Subjects with Serious Adverse Event, n	Subjects with Serious Adverse Events, %	Serious Adverse Events, n	Serious Adverse Events per Subject-Year
Angina	1	0.6%	1	0.019
Arrhythmia	2	1.3%	2	0.038
Cardiac Arrest	1	0.6%	1	0.019
Infection	6	3.9%	6	0.115
Influenza	1	0.6%	1	0.019
MRI Incidental Finding	1	0.6%	1	0.019
Medication Related	2	1.3%	2	0.038
Medication Related CHF	1	0.6%	1	0.019
Musculoskeletal	2	1.3%	2	0.038
Musculoskeletal Injury	1	0.6%	1	0.019
Orthopedic Injury	1	0.6%	1	0.019
Peripheral Vascular Disease	1	0.6%	1	0.019
Psychiatric	1	0.6%	1	0.019
Renal	2	1.3%	3	0.057
Respiratory	2	1.3%	2	0.038
Stroke	1	0.6%	1	0.019
Syncope/Pre-Syncope	1	0.6%	1	0.019
Worsening CHF	9	5.8%	12	0.229
Total	23	14.9%	40	0.765

Number of Enrolled Subjects = 154, Number of Subject-Years since Enrollment = 52.3

Table 17: Summary of Non-serious Adverse Events

Adverse Event	Subjects with Adverse Event, n	Subjects with Adverse Events, %	Adverse Events, n	Adverse Events per Subject-Year
Arrhythmia	5	3.2%	5	0.096
Device Related	2	1.3%	2	0.038
Gastrointestinal	1	0.6%	1	0.019
Infection	3	1.9%	3	0.057
MRI Incidental Finding	24	15.6%	25	0.478
MRI Procedure Discomfort	6	3.9%	6	0.115
Medication Induced CHF	1	0.6%	1	0.019
Medication Related	3	1.9%	3	0.057
Musculoskeletal	6	3.9%	6	0.115
Neurological	1	0.6%	1	0.019
Other	9	5.8%	12	0.229
Other - MRI Procedure	1	0.6%	1	0.019
Other - Medical	1	0.6%	1	0.019
Peripheral Vascular Disease	1	0.6%	1	0.019
Pocket Pain/Discomfort	4	2.6%	5	0.096
Psychiatric	1	0.6%	1	0.019
Renal	1	0.6%	1	0.019
Respiratory	3	1.9%	3	0.057
Syncope/Pre-Syncope	1	0.6%	1	0.019
Worsening CHF	5	3.2%	5	0.096
Total	57	37.0%	84	1.625

Number of Enrolled Subjects = 154, Number of Subject-Years since Enrollment = 52.3

ProMRI® System

Technical Manual



Manufactured by:

BIOTRONIK SE & Co. KG
Woermannkehe 1
12359 Berlin Germany

M4188-C 08/16
© 2016 BIOTRONIK, Inc. All rights reserved.
MN062r4 8/29/16

www.biotronikusa.com/promri

BIOTRONIK, Inc.
6024 Jean Road
Lake Oswego, OR 97035-
5369

(800) 547-0394 (24-hour)
(800) 291-0470 (fax)
www.biotronik.com

 **BIOTRONIK**
excellence for life