


Acticor 7 HF-T/HF-T QP Product Highlights

NEW



Therapy Highlights

CRT AutoAdapt

Continuous AV delay optimization and automatic pacing chamber selection

AV AutoAdapt

Continuous AV delay optimization automatically adjusting to changes in patient condition

Up to 20 LV Pacing Vectors

Industry-leading choice of vectors to optimize CRT response and address complications

Auto LV VectorOpt

Automatic pacing threshold testing with pre-selection options and direct programming

AV Opt

Facilitated AV delay optimization via the programmer

Device Hardware Highlights

BIOshape

Reduced skin pressure and increased patient comfort

Up to 9 years* longevity

Fewer box replacements leading to reduced risks, complications and costs

*Parameters for HF-T QP: 60 ppm; RA 15%, RV/LV 100% pacing, RA/RV/LV @ 2.5 V/0.4 ms; 500 ohm, 2 max. energy shocks/year. Home Monitoring ON; diagnostics ON.

DX Technology with DF-4

Atrial diagnostics without an atrial lead for reduced risks, complications and costs

ProMRI Highlights

3T Full-Body MRI Scanning with MRI AutoDetect

High image quality and reduced scanning time now available for all parts of the body

BIOTRONIK Home Monitoring Highlights

QuickCheck

Fast access to patient data (3-4 min) for more efficient treatment and greater peace of mind for patients

More Life. Made Simpler

 **BIOTRONIK**
excellence for life