

# Acticor 7 VR-T, VR-T DX and DR-T Product Highlights

NEW



BIOshape

Ultraslim

## Device Hardware Highlights

### BIOshape

Reduced skin pressure and increased patient comfort

### Up to 15 years\* longevity (VR-ICD)

Fewer box replacements leading to reduced risks, complications and costs

\*Parameters for VR-T: 40 ppm, RV 15% pacing, RV 2.5V/0.4 ms; 500 ohm; Includes 2 max. energy shocks/year; Home Monitoring ON; diagnostics ON.

## Therapy Highlights

### Auto Far-Field Protection (FFP)

Automatic shortening of the Far-Field Protection interval after VS to avoid functional atrial undersensing

### SMART Algorithm Enhancements

Improved differentiation between VT with retrograde conduction and 1:1 SVTs (AT, AVNRT, AVRT) or atrial flutter

### ProMRI Highlights

#### 3T Full-Body MRI Scanning with MRI AutoDetect

High image quality and reduced scanning time now available for all parts of the body

## Diagnostics Highlights

### Extended Lead Measurement (XLM)

Extended range of impedance measurements for a more sophisticated lead performance assessment

### DX Technology with DF-4

Atrial diagnostics without an atrial lead for reduced risks, complications and costs

## BIOTRONIK Home Monitoring Highlights

### QuickCheck

Fast access to patient data (3-4 min) for more efficient treatment and greater peace of mind for patients

More Life. Made Simpler

 **BIOTRONIK**  
excellence for life